

7MINUTES PLASMA STERILIZER

PLASMAPP

(저온플라즈마 7분멸균기)

Low_Temperature_PLASMA_Sterilizer

50 Countries – 1500 References



FPS-15s Plus(14L)



MINI(7L)

Plasmapp Sterilization Product Catalog

Fast Low-Temperature Plasma Sterilizer – 55°C



STERLINK®

Chamber : 14 Liter
3 Cycles Mode - Pouch : 7 min
- Pouch+ : 14 min
- Full Load : 36 min
Class IIb, CE Certified



STERLINK® Mini

Chamber : 7 Liter
2 Cycles Mode - Pouch : 7 min
- Full Load : 18 min
Class IIb, CE Certified

Specifications



Model	Dimensions [mm]			Weight [kg]	Chamber [liter]	Cassette (Sterilization cycle time)		
	Width	Depth	Height			STERPACK®	STERPACK® Plus	STERLOAD®
STERLINK®	433	614	437	67	14	○ (7 min)	○ (14 min)	○ (36 min)
STERLINK® Mini	270	450	330	18	7	○ (7 min)	X	○ (18 min)

PLASMAAPP

(LOW_TEMPERATURE_STERILIZER)

(세계최초 직분사방식 특허기술 적용,
안전멸균율 : $SAL10^{-6}$ 사전검증)

빠른 멸균
저온 멸균
깨끗한 멸균

안전 멸균율(SAL:99.9999%)

저온 7분 멸균
컴팩트한 디자인
원터치 운전
실시간 모니터링

적용분야 & 잠재고객

실험실, 연구소, 대학교, 바이오메디컬 산업분야, 제약회사,
종합병원(안과, 치과 제외)카지노, 수의과, 동물병원 등

1. FPS-15s Plus 사양

분 류	내 용			단 위
치 수	433 (W) x 614 (D) x 437 (H)			mm
챔버 용량	14			L
안전장치	물 리 적	도어센서 및 진공에 의한 열림 방지		
	전 기 적	과전류 방지용 정격 퓨즈 삽입		
공정변수	최대온도	60 미만		℃
	압력범위	0 – 760		torr
	과산화수소 주입 직전 압력	3.0 이하		torr
	확산압력	20 이상, 80 이하		torr
	공정확인	Chemical indicator (CI)	Class 4 (Compliance: ISO 11140-1)	
진 공	최저압력	3 이하		torr
	압력변화율	60 초 동안 6 torr 이하의 압력변화를 확인		
준비공정성능 (Smart Ready™ Process)	바코드	자동 바코드 인식 및 MODE 선택		
	건조성능	최저압력	3.0 이하	torr
		최저 상대습도	5 이하	%
멸 균	멸균성능	Biological indicator (BI)	SCBI (Compliance: ISO 11138)	
		LUMEN (Single-channel SUS)	Dia: 0.7, Length: 500	mm
		LUMEN (Single-channel SUS)	Dia: 2.0, Length: 1,500	mm
		LUMEN (Single-channel PTFE)	Dia: 1.0, Length: 2,000	mm
	멸균시간	POUCH	4	min
		POUCH PLUS	8	
		CHAMBER	15	
	Full cycle 소요시간 (기본 사양)	POUCH	7 ± 1 이내	
		POUCH PLUS	14 ± 1 이내	
		CHAMBER	36 ± 1 이내	
	Full cycle 소요시간 (최대 적재)	POUCH	10 ± 1 이내	
		POUCH PLUS	18 ± 1 이내	
멸균제	과산화수소수 (H ₂ O ₂ 58 – 59.5%)			
정화방식	플라즈마 방전 처리 (Low frequency DBD)			
Display	치 수	153.1 (W) x 85.6 (H)		mm
	해 상 도	800 x 480		pixels
포트	IEC-320 C14, USB type B, RS-232, RJ-45			
전기적 특성	정격전원	220-240 (10 A)		VAC
	정격주파수	50/60		Hz

2. STERLINK MINI 사양

분 류	내 용			단 위
치 수	본체: 270 (W) x 450 (D) x 330 (H) 펌프 유닛: 270 (W) x 450 (D) x 330 (H)			mm
챔버 용량	7			L
안전장치	물 리 적	도어센서 및 진공에 의한 열림 방지		
	전 기 적	과전류 방지용 정격 퓨즈 삽입		
공정변수	최대온도	60 이하		℃
	압력범위	0 - 760		torr
	과산화수소 주입 직전 압력	3.0 이하		torr
	확산압력	20 이상, 80 이하		torr
	공정확인	Chemical indicator (CI)	Class 4 (Compliance: ISO 11140-1)	
진 공	최저압력	3.0 이하		torr
	압력변화율	0.1 이하		torr·s ⁻¹
멸균공정 (Sterilization) ,60Hz	멸균성능	LUMEN (Single-channel Stainless Steel)	Dia: 0.7, Length: 500 Dia: 2.0, Length: 1,500	mm
		LUMEN (Single-channel PTFE)	Dia: 1.0, Length: 2,000	mm
	멸균시간	POUCH	4	분
		CHAMBER	8	
순환시간 (Cycle time) ,60Hz	반주기 공정 (Half cycle)	POUCH	5	분
		CHAMBER	13	
	전주기 공정 (Full cycle)	POUCH	7	분
		CHAMBER	18	
멸균제	과산화수소수 (H ₂ O ₂ 58 - 59.5%)			
정화방식	플라즈마 방전 처리 (Low frequency)			
Display	치 수	153.1 (W) x 85.6 (H)		mm
	해 상 도	800 x 480		pixels
포트	IEC-320 C14, RJ-45, USB type A, USB type B, cable connetor			

YOU TUBE :

1. What is the Plasmapp ? :

<https://youtu.be/NpcH8QGYUc4>

2. Plasmapp_Low temperature Plasma Sterilizer_FPS-15s :

<https://youtu.be/OfAMkMkl7Xk>

(New)Plasmapp_Low temperature Plasma Sterilizer_FPS-15s :

<https://youtu.be/-mZC3OVOR9A>

3. Sterlink :

<https://youtu.be/fk3tF7WTt9I>

4. Introduction Plasmapp :

<https://youtu.be/sQeYH0u1w40>

TEST-INSPECTION REPORT

Mask Sterilization Test

Manufacturer Name: Plasmapp Co., Ltd.
Representative: You Bong, LIM
Location: 83 Jukdong-ro, Yuseong-gu, Daejeon 34127, Republic of Korea

Product Name: Low temperature plasma sterilizer
Brand Name: STERLINK®
Model Name: STERLINK MINI
Serial Number: M07BUG042A

Test-Inspection Item: Mask Sterilization Test
Testing Laboratory: Plasmapp Research Institute
Location: 327 Dongbu-daero, Osan-si, Gyeonggi-do 18151, Republic of Korea (South Korea)

Decision: Pass

Tested until: 27 Apr. 2020 **Issued Date:** 29 Apr. 2020
Tested by: Hyun Jeong, JEON **Approved by:** Seung Hun, LEE
Process Inspector General Director of R&D



Plasmapp Research Institute



Mask Sterilization Test

1. Test schedule

1.1 Date of test beginning: 14 Apr. 2020
1.2 Date of test completion: 27 Apr. 2020

2. Test article

Low temperature plasma sterilizer (STERLINK® MINI, S/N: M07BUG042A)

3. Test guideline

3.1 The tests were performed in accordance with the standard of ISO 14937:2009.

3.2 Information of testing materials

3.2.1 Masks

Type	Product	Model	Type	Manufacturer	Expiration date
KF94	All guard basic dust mask	-	KF94	All guard	Mar. 2023
N95	Particulate Respirator N95	81105	N95	3M	Dec. 2024

3.2.2 Bacteria spore suspension

Item	Details
Model	Spore suspension
Manufacturer	Mesal abs
Lot Number	SSS-847
Expiration date	25 Oct. 2021
Comments	1.3 × 10 ⁶ <i>Geobacillus stearothermophilus</i> spores (ATCC Cell Line 7953) / 0.1 ml complying with ISO 11138

*Figures of the masks and inoculated position were referred to the Appendix 1.

1/7

3.3 Test methods

3.3.1 Half cycle sterilization test

- (1) The *Geobacillus stearothermophilus* spore suspension was inoculated to the specified area of each mask as shown in the Appendix 1 and dried. Before the inoculation, the spore suspension for inoculation was diluted by serial 10-fold dilution method and they were inoculated in tryptic soy agar as triplicate plates for each dilution to confirm the number of the spores. The inoculated plates were incubated at 58°C for 48 hours. The prepared spore suspension for the inoculation should have greater than 1.0 × 10⁶ colony forming unit (CFU).
- (2) Each mask with spores was inserted in Tyvek® pouch and sealed. **For the worst-case condition, the masks were positioned as far as possible from the sterilant injection nozzle.**
- (3) The prepared pouch was processed with half cycle sterilization of chamber mode.
- (4) For the control, inoculated mask which was left on the room temperature without sterilization process was used.
- (5) After sterilization cycle, the inoculated area on the mask was cut by a flame-sterilized scissor and placed in the tryptic soy broth of 5 ml.
- (6) All test and control sample were incubated at 58°C for 7 days. After incubation, all samples were mixed by a vortex mixer and the turbidity change of media was inspected. The results were reported as positive (growth) or negative (no growth).
- (7) The mask sterilization test except control was performed in three consecutive half sterilization cycles.

2/7

4. Test results**

4.1 The number of the spores for inoculation

	Number of recovered spores (10 ⁶ CFU)
#1	1.17 ± 0.08

4.2 Results of the mask sterilization test

Mask type	Condition	Number of positive/Number of tested			Sterilization Results
		#1	#2	#3	
KF94	Control		1/1		-
	Test	0/1	0/1	0/1	Pass
N95	Control		1/1		-
	Test	0/1	0/1	0/1	Pass

5. Conclusions

- (1) The number of spores used for inoculation was above 1.0 × 10⁶ CFU.
- (2) After three consecutive mask sterilization tests, all the incubation results were negative except control.
- (3) According to the test results, mask sterilization test is completely successful.

**The related figures were referred to the Appendix 2. The time evolution of pressures inside the chamber were described in the Appendix 3, as well.

3/7

Appendix 1

1. Masks




Figure 1.1 (a, c) The outside and (b, d) the inside for KF94 and N95, respectively. The area of inoculation was indicated as the square filled with diagonal line.

4/7

2. Half cycle of chamber mode for the sterilization test




Figure 1.2 (a) The position of the inoculated masks sealed by Tyvek® pouch. The figure of (b, d) before and (c, e) after the sterilization for KF94 and N95, respectively.

Appendix 2

1. The number of the spores for the mask inoculation




Figure 2.1 Incubation results of 10⁴ dilution samples of the spore suspension for inoculation.

2. Results of the mask sterilization test




Figure 2.2 The incubation results for mask sterilization test, (a, c) The control and (b, d) three consecutive mask sterilization tests for KF94 and N95 type masks, respectively. The turbidity can be determined by black lines behind the tubes.

Appendix 3

1. Pressure curves of the mask sterilization test



Figure 3.1 (a) The whole plot of pressure curve during the half cycle and (b) magnified plot of diffusion phase for each consecutive trial



BEYOND ASIAN HUB, TOWARDS GLOBAL WORLD

TEST REPORT

98, Gyojuukwon-ro, Gwacheon-si, Gyeonggi-do, 13810, Korea TEL 82-2-2164-9011 FAX 82-2-2634-1008
Report No. : MSK-2020-00064 Receipt Date : 2020.03.18
Representative : Youdong Lim Test Completion Date : 2020.04.20
Company name : Plasmapp Co., Ltd.
Address : 1F, 3F, #301, 83 Jukdong-ro, Yuseong-gu, Daejeon 34127
Sample name : Sterilizer, plasma, low temperature (Serial No. P155UB0315A)

Test Results

TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
Virucidal test (Human Influenza A (H1N1))	Log reduction	PPS-15s Plus	33.38	By the client (Annex#1)

** Test virus

- Human Influenza A (H1N1) (Strain: KJMC-76), KBPV-VR-76

*** Host cells

- MDCK (NB1-2), ATCC CCL-34

**** Log reduction (LR)

- LU = LU - LT

- LU = The titer recovered from the virus control specimen (untreated)

- LT = The titer recovered from virucidal test specimen (treated)

***** Interpretation of the results

- log reduction 1) 1.38 %

- log reduction 2) 2.00 %

- log reduction 3) 3.68 %

- log reduction 4) 4.69 %

- log reduction 5) 5.69 %

- Test period : Mar.18,2020 ~ Apr.20,2020

Annex#1 : By the client - The test methods

Annex#2 : Picture of appearance or shape

Attachment : Test Report

- Next Page -

Hyun Ho Kim
Prepared by Hyun Ho Kim
Tel : 02-2662-3998

Yin-Sik Cho
Reviewed by Jin-Sik Cho
Technical Manager
Tel : 02-2661-0955 (3~4)

Korea Testing & Research Institute

President *Kwon Oh-jung*



QR Code for forgery

1 of Total 2 Page(s)

KTR KOREA TESTING & RESEARCH INSTITUTE



BEYOND ASIAN HUB, TOWARDS GLOBAL WORLD

TEST REPORT

98, Gyojuukwon-ro, Gwacheon-si, Gyeonggi-do, 13810, Korea TEL 82-2-2164-9011 FAX 82-2-2634-1008
Report No. : MSK-2020-00064 Receipt Date : 2020.03.18
Representative : Youdong Lim Test Completion Date : 2020.04.20
Company name : Plasmapp Co., Ltd.
Address : 1F, 3F, #301, 83 Jukdong-ro, Yuseong-gu, Daejeon 34127
Sample name : Sterilizer, plasma, low temperature (Serial No. P155UB0315A)

Test Results

TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
- Usage of Report : QUALITY CONTROL				

NOTE : 1. The test results of this test report are only limited in to the samples and sample names provided by the client and do not guarantee the quality of all products of the client. You can check website (www.ktr.or.kr) or QR code to verify the authenticity of the certificate.
2. This test report shall be used only within the purpose of its defined usage and shall not be used to public relation, advertisement and lawsuit.
3. This test report is only valid when printed on KTR original report paper with hologram and when re-issued by KTR. The copy and the electronic file of the test report are only for reference.

Hyun Ho Kim
Prepared by Hyun Ho Kim
Tel : 02-2662-3998

Yin-Sik Cho
Reviewed by Jin-Sik Cho
Technical Manager
Tel : 02-2661-0955 (3~4)

Korea Testing & Research Institute

President *Kwon Oh-jung*



QR Code for forgery

2 of Total 2 Page(s)

KTR KOREA TESTING & RESEARCH INSTITUTE

Annex#1 : By the client – The test methods

1. Test name

Virucidal test

2. Test article

Product name: Low temperature plasma sterilizer
Brand name: STERLINK®
Model name: FPS-15s Plus (FPS-15s)
Consumable (Pouch with sterilant): STERPACK®

3. Virus used for test

Influenza A (H1N1)

4. Test methods

- 1) The virus was inoculated to the sterilized glass petri dish (5cm diameter) and dried.
- 2) The sterilized glass petri dish with virus was inserted in the STERPACK® and sealed.
- 3) The prepared STERPACK® was installed to the loading block of FPS-15s Plus and processed with half cycle sterilization of pouch mode.
- 4) After the sterilization, the virus inoculated on the sterilized glass petri dish was recovered and inoculated to the host cell to confirm the virus titer.
- 5) For the control, the sterilized glass petri dish with virus was inserted in the STERPACK®, sealed, and left on the room temperature without sterilization process during the sterilization time.
- 6) The log reduction was calculated by comparing the test results of the test and control sample. According to the test results, the virucidal performance of the FPS-15 Plus was confirmed.

Annex#2 : Picture of appearance or shape



KTR

Test Report

MSK-2020-000694

Low temperature plasma sterilizer

Virucidal test

President of KTR

Test outline

Study title : Virucidal test

Study number : MSK-2020-000694

Sponsor

Name : Plasmapp Co., Ltd.
Address : 83 Jukdong-ro, Yuseong-gu, Daejeon 34127, Republic of Korea
Representative : You Bong, LM

Test facility

Name : Korea Testing And Research Institute
Address : 98, Gyooyukwon-ro, Gwacheon-si, Gyeonggi-do, Korea

Hyun Ho Kim 2020.04.20.
Date
Hyun Ho Kim, Ph.D.
Study Director
Medical device-Bio Research Institute, KTR
Cho Jin-Sik 2020.04.20.
Date
Cho, Jin-Sik, M.S.
Technical Manager
Medical device-Bio Research Institute, KTR

This report is presented for the article that was submitted by the sponsor.

April 20, 2020

CONTENTS

1. Summary	1
1.1. Test Schedule	1
1.2. Picture of the Test Article	1
2. Equipment & Materials	2
2.1. Test Equipment	2
2.2. Test Article	2
2.3. Test Materials	2
3. Test Method	4
3.1. Host Cell Culture	4
3.2. Viral Culture	4
3.3. Viral Film Formation	4
3.4. Virucidal Test	4
3.5. Validation Test	5
3.6. Virus Recovery	5
3.7. Calculation of Results	6
3.8. Test Requirements	6
4. Results	8
4.1. Cytotoxicity of the Test Solution	8
4.2. Results of Virucidal Test	8
4.3. Results of Validation Test	9
5. Conclusion	10
6. References	10

Tables

Table 1. Results of the Cytotoxicity Test	8
Table 2. Results of Influenza A (H1N1)	8
Table 3. Results of Calculation	8
Table 4. Results of Log Reduction	9
Table 5. Results of the Sub-cytotoxicity and Neutralization Test	9

1. Summary

This test was performed to assess the virucidal activity of the test article provided by the sponsor [Low temperature plasma sterilizer (Model: FPS-15s Plus)] according to By the client, The Human Influenza A (H1N1) requested by the sponsor were used as test viruses in the test. The virus suspension was dried on an inanimate, nonporous surface. After placing the viral film inside the test Pouch requested by the sponsor and operating the article, and then the eluates are assayed for infectivity. The activity of the virus was quantified by tissue culture 50 % infectious dose assay (log₁₀TCID₅₀) after infecting the host cells with the virus.

Under the test conditions, the log reduction of the Influenza A (H1N1) in the test article [Low temperature plasma sterilizer] was found to be >3.38 after operation of the test article.

1.1. Test Schedule

Total test period March 18, 2020 to April 20 2020

1.2. Picture of the Test Article



MSK-2020-000694

2. Equipments & materials

2.1. Test Equipments

Biosafety cabinet class II	(Cryste, Republic of Korea)
Autoclave	(Coretech, Republic of Korea)
Dry oven	(Jisico, Republic of Korea)
Water bath	(Fisher, USA)
CO ₂ Incubator	(Thermo, USA)
Centrifuge	(Hani, Republic of Korea)
Deep freezer	(Thermo, USA)
Microscopy	(Leica, USA)
pH meter	(Thermo Orion, USA)
Vortex mixer	(Thermolyne, USA)
Haemocytometer	(Marienfeld, Germany)
75-cm ² flask	(Falcon, USA)
150-cm ² flask	(Falcon, USA)
Test tube	(Falcon, USA)
Serological pipette	(Falcon, USA)
96 well plate	(Falcon, USA)
Cell scraper	(Falcon, USA)
0.45 μm filter	(Advantec, USA)

2.2. Test Article

Low temperature plasma sterilizer, Model: FPS-15s Plus, Serial No.: P15BU015A
(provided by the client)

Pouch, Model: STERPACK®, Serial No.: SP20C002
(provided by the client)

Glass petri dish, Model: -, Lot No.: -
(provided by the client)

2.3. Test Materials

2.3.1. Virus
Human Influenza A (H1N1) (KBPV-VR-76)

2.3.2. Host Cells
MDCK (NBL-2) (ATCC CCL-34)

2.3.3. Media and Reagents

Minimal Essential Media (MEM)	(GIBCO, USA)
Phosphate Buffered Saline (PBS)	(GIBCO, USA)
Fetal Bovine Serum (FBS)	(GIBCO, USA)
Bovine Serum Albumin (BSA)	(Sigma, USA)
Penicillin-Streptomycin (PS)	(GIBCO, USA)
Trypsin-EDTA	(GIBCO, USA)
TPCK-treated Trypsin	(Sigma, USA)
Crystal Violet	(Sigma, USA)
Methyl Alcohol	(Sigma, USA)

2.3.4. Cell Culture Medium
MEM + 10 % FBS

2.3.5. Virus culture medium
MEM + 0.3 % BSA + 1 μg/mL TPCK-treated trypsin

3. Test method

3.1. Host Cell Culture

MDCK cells were cultured by using a cell culture medium under the conditions of 5 % CO₂ and (36 ± 2) °C.

3.2. Viral Culture

3.2.1. Influenza A (H1N1) Culture
MDCK cells monolayer-cultured in a 75-cm² flask were washed with PBS, and then infected with 3 mL of Influenza A (H1N1) diluent. After infection with the virus as described in section 3.2.2, the virus was cultured for 3 to 5 days until more than 90 % of host cells showed a cytopathic effect by adding virus culture medium. Subsequently, the culture supernatant was centrifuged at 2000 rpm for 10 minutes and the supernatant was filtered with a 0.45 μm filter to isolate the virus.

3.2.2. Virus treatment conditions

Virus	CO ₂ concentration	Temperature	Infection time
Influenza A (H1N1)	5 %	(33 ± 2) °C	(60 ~ 90) minutes

3.3. Viral Film Formation

After adding 0.2 mL of the test virus to the sterilized glass petri dish, the viral film was formed by drying for 30 minutes in the ambient conditions of the laminar flow BSC with the petri dish cover removed.

3.4. Virucidal Test

3.4.1. Cytotoxicity Test
As in Section 3.3., 0.2 mL of the virus culture medium was placed on the inside bottom surface of a glass petri dish to form a film. The film was placed inside the test pouch requested by the client and sealed, and the test article was operated in a HALF CYCLE(operation time: about 6 minutes). After overlaying 2 mL of the virus culture medium, the film was resuspended by scraping the surface of the glass petri dish using a cell scraper(the considered the 10¹ dilution, hereinafter the same). And then the mixture was serially diluted with the virus culture medium by a factor of 10. The diluted solutions were treated with the host cells at (36 ± 1) °C for (30 ~ 60) minutes, and the cytotoxicity effect was observed by using a microscope.

3.4.2. Cell Culture Control

The MDCK cells were treated with the virus culture medium, which was used as the cell culture control.

3.4.3. Virus Susceptibility Control

The test virus was mixed at 1:9 ratio in virus culture medium, and diluted stepwise by a factor of 10 in the virus culture medium. Each diluent was infected into the host cells to measure the viral titer of the test virus(initial).

After forming the viral film as in Section 3.3., it was sealed in a test Pouch requested by the sponsor and allowed to stand at (22 ± 2) °C for 6 minutes. Thereafter, 2 mL of the virus culture medium was added, and then the virus film was recovered by scraping the surface of the glass petri dish using a cell scraper, and then diluted 10 times in the virus culture medium. Each diluent was infected into the host cells to measure the viral titer of the control group(After 6 minutes).

3.4.4. Virus Test

After forming the viral film as in Section 3.3., it was sealed in a test Pouch requested by the sponsor, and the test article was operated in a HALF CYCLE(operation time: about 6 minutes). Subsequently, the viral titer of the test group(After operation) was measured in the same manner as in Section 3.4.3.(the viral titer of the control group).

3.5. Validation Test

The solution recovered as in the cytotoxicity test(Section 3.4.1.) was used as the test solution, and used for the validation test (below).

3.5.1. Sub-cytotoxicity Test

The test solution for each concentration not exhibiting cytotoxicity was treated with the host cells for 30 minutes, and then the cell monolayer was washed with PBS. Subsequently, the cells were infected with the test virus dilution as described in section 3.2.2, and cultured as described in section 3.6.1..

3.5.2. Neutralization control

The neutralized test solution was mixed with the same amount of the test virus dilution, and the mixture was kept for (10 ~ 20) minutes for reaction. Subsequently, the cells were infected with the mixture as described in section 3.2.2., and cultured as described in section 3.6.1..

3.6. Virus Recovery

After the reaction, each solution was serially diluted 10-fold, and 100 μL of each diluted solution treated in 8 wells (8 replicates) of the cells prepared in advance in 96 well plates. Then, the cells were incubated as described in section 3.2.2., and cultured as described in section 3.6.1..

3.6.1. Culture Conditions

Virus	CO ₂ concentration	Temperature	Culture time
Influenza A (H1N1)	5 %	(33 ± 2) °C	(3 ~ 5) days

3.6.2. Staining Conditions

5. Conclusions

- (1) The results of bacteriostasis tests were shown as all positive except negative controls.
- (2) According to the bacteriostasis test results, the growth of *Geobacillus stearothermophilus* spore inoculated on the mask were not inhibited by three consecutive full sterilization cycle of chamber mode.

Appendix 1

1. Masks

Figure 1.1 (a, c) The outside and (b, d) the inside for KF94 and N95, respectively. The area of inoculation was indicated as the square filled with diagonal line.

2. Chamber mode for the bacteriostasis test



Figure 1.2 (a) The position of the masks sealed by Tyvek® pouch. The figure of (b, d) before and (c, e) after the sterilization for KF94 and N95, respectively.

Appendix 2

1. The number of the spores in spore suspension prepared for inoculation.

Figure 2.1 The incubation results of 0.1 ml spore suspension prepared for bacteriostasis test. They were incubated at 58 °C for 24 hours.

Appendix 3

1. Results of the bacteriostasis test




Figure 3.1 The incubation results for bacteriostasis test. (a) The positive control, (b) The negative control and (c) test sample of the KF94 type mask. (d) The negative control and (e) test sample of the N95 type mask. The turbidity can be determined by black lines behind the tubes.

Appendix 4

1. Pressure curves of the bacteriostasis test

Figure 4.1 (a) The whole plot of pressure curve during the full cycle and (b) magnified plot of diffusion phase for each consecutive trial














US Veterinary Clinics



Spain Ophthalmology Clinics



France Dental Clinics



India Ophthalmology Clinics



Hong Kong Dental Clinics



Germany Dental Clinics




USA Partner company



USA Vet Hospital – Malibu, CA



USA Vet Hospital – Sierra, NV



USA Vet Hospital – Palisade, CA





STERLINK®

7 Minutes Low-Temperature Plasma Sterilizer



 **Plasmapp**

STERLINK[®]



Fast 7 min Sterilization



Compact Design



One-touch Operation



Real-time Monitoring



STERLINK® FPS-15s Plus 14L | STERLINK® MINI 7L


A low-temperature plasma sterilizer Sterilization of a wide range of medical devices

The world's first direct injection pouch-type
patented sterilization technology has been applied.
SAL of 10^{-6} is validated to ensure sterilization stability of medical devices.

(SAL: Sterility Assurance Level, 99.9999%)



Smart Sterilization



Sterilant Cassettes, Consumables and Accessories

※ Sterilant cassettes



STERPACK®
· Size : 135×280mm
· 50ea / Box



STERPACK® Plus
· Size : 240×410mm
· 30ea / Box



STERLOAD®
· 30ea / Box



STERLOAD® MINI
· 30ea / Box

※ Consumables



Tyvek® 100
· Size : 100×400mm
· 120ea / Box



Tyvek® 200
· Size : 200×400mm
· 90ea / Box



Tyvek® 300
· Size : 300×400mm
· 60ea / Box



Sterilization wrap
· Size : 60×60cm
· 50ea / Box

Accessories



STERSEAL®
· Size : 505×255×145mm
· Weight : 12kg
· Temperature : 135℃ / 150℃



Cart
· Size : 483×660×603mm
· Weight : 37kg



Label printer
· Size : 120×102×146mm
· Weight : 0.5kg



Label sticker roll
· Width : 52 × 110mm
· Label counts : 49ea

※ Accessories



CI tape
· Width : 20mm
· Length : 50m



CI strip
· Size : 18×105mm
· 250ea / Box



BI
· Time : 30 min
· 50ea / Box

Tray



STERPACK® Tray
· Size : 195×80×30mm



STERPACK® Plus Tray
· Size : 260×160×50mm

※ Single Use Only

※ Tyvek® is a Dupont™ registered trademark.



STERPACK®
STERLOAD®



The world's first direct injection pouch with patented sterilization technology



Vacuum sealing pouch enables to check immediate sterilization conditions



Using safe disposable sterilant cassettes

STERPACK® Specification

Cassette	Mode	Sterilant	Cycle time	Capacity
STERPACK®	Pouch Mode	H ₂ O ₂ 58%–59.5% (0.1 cc/cell)	7min	1L
STERPACK® Plus	Pouch Plus Mode	H ₂ O ₂ 58%–59.5% (0.3 cc/cell)	14min	4L

STERLOAD® Specification

Cassette	Mode	Sterilant	Cycle time	Capacity
STERLOAD®	Chamber Mode	H ₂ O ₂ 58%–59.5% (0.9 cc/cell) (0.7cc/cell)	36min	14L
STERLOAD® MINI			18min	7L

※ STERPACK® Plus, STERLOAD® : Compatible with STERLINK® FPS-15s Plus only

※ STERLOAD® MINI : Compatible with STERLINK® MINI only

Using Tyvek® pouch



Using in chamber mode




Sterilization Performance

STERLINK® verifies sterilization performance through the following lumen tests

Lumen claims for STERLINK®


- 0,7ø x 500 mm Stainless
- 2,0ø x 1500 mm Stainless
- 1,0ø x 2000 mm PTFE



Instrument Tracking System (ITS™)

Equipment tracking using barcode

- Real-time monitoring to check the STERLINK®
- Provide the latest software updates remotely



Comparison Chart

Sterilizer	Temperature	Cycle Time	Sterilant	Feature
Autoclave	Up to 134°C	60min + 1 hour cooldown	Hot steam	• Cloth sterilization • Long cycle time • Risk of burn • Damage of medical instruments
E.O. Gas	Up to 60°C	Over 120min + 12 hour ventilation	E.O. gas	• Highly dangerous toxic gas • Long cycle time
Plasma Sterilizer	Up to 60°C	70min	H ₂ O ₂	• High cost • Large volume • Additional ventilation • Very low efficiency
STERLINK®	Up to 60°C	Pouch Mode : 7min Pouch Plus Mode : 14min Chamber Mode : 36min	H ₂ O ₂ direct injection	• Fast sterilization cycle • Economic cost • Compact size • Ergonomic design • Easy maintenance • ITS™ system • Eco-friendly • Safeness of heat sensitive instrument

STERLINK®

Specification



Model	STERLINK® FPS-15s Plus	STERLINK® MINI
Size (W×D×H)	433 × 614 × 437 mm	275 × 440 × 330 mm
Chamber (W×D×H)	264 × 410 × 125 mm (14L)	206 × 346 × 113 mm (7L)
Diagonal length (Chamber)	47cm	36cm
Vacuum Pump	Pump built-in type	Pump stand-alone type
Weight	67kg	20kg (Pump module : 21kg)
Cycle time	STERPACK® : 7min STERPACK® Plus : 14min STERLOAD® : 36min	STERPACK® : 7min STERLOAD® MINI : 18min

Applications & potential customers

Laboratory, Research Center, University, Bio
Medical Industries, Pharmaceutical Company,
General Hospital (except
Ophthalmology/Dental clinic), Casino,
Veterinary medicine, Animal Hospital etc.,

Contact : plasmapp.kr@gmail.com

